

Anti-Phospho-P53 (Ser15) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P04637

Item Number: CM0001774



Introduction

KO-validated rabbit monoclonal antibody against human phospho-p53 (Ser15). Suitable for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse p53. Ideal for cancer and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Phospho-p53-S15 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S15
Species	Human
Gene ID (Human)	7157
Gene ID	7157
Immunogen	A synthetic phosphorylated peptide around S15 of human p53
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04637
Protein Weight	44 kDa
Shipping	Ice bag