


Anti-Phospho-P53 (T55) Rabbit Monoclonal Antibody For Wb And Elisa - P04637

Item Number: CM0001095



Introduction

Phospho-p53-T55 Rabbit mAb (Cat. AP0986) validated for WB and ELISA. Specifically detects human p53 phosphorylated at Thr55, enabling studies of p53 activation and tumor suppression mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-T55 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-T55
Species	Human
GeneID	7157
Immunogen	A synthetic phosphorylated peptide around T55 of human p53 (P04637)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag