


Anti-Phospho-P53 (T81) Rabbit Polyclonal Antibody For Wb, Elisa - P04637

Item Number: CM0001836



Introduction

Rabbit polyclonal antibody targeting phosphorylated p53 at Thr81, validated for WB and ELISA in human and rat. Ideal for studying apoptosis, cell cycle regulation, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-T81 Rabbit pAb
Remarks/Alias	P53, BCC7, LFS1, BMF55, TRP53, Phospho-p53-T81
Species	Human
GeneID (human)	7157
GeneID	7157
Immunogen	A synthetic phosphorylated peptide around T81 of human TP53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag