

Anti-Phospho-P70 S6 Kinase 1-S371 Rabbit Polyclonal Antibody For Wb And Elisa - P23443

Item Number: CM0000752



Introduction

Rabbit polyclonal antibody detecting p70 S6 kinase 1 phosphorylated at Ser371. Validated for WB and ELISA; reacts with human and mouse RPS6KB1. Supports mTOR signaling and translation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-S371 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-S371
Species	Human
GeneID (Human)	6198
GeneID	6198
Immunogen	A synthetic phosphorylated peptide around S371 of human RPS6KB1 (NP_001258989.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag