

Anti-Phospho-P70 S6 Kinase 1 (Ser371) Rabbit Pab For Wb And Elisa - P23443

Item Number: CM0000427



Introduction

High-affinity rabbit polyclonal antibody specifically recognizing human p70 S6 Kinase 1 phosphorylated at Ser371. Suitable for Western blot (WB) and ELISA; cross-reacts with mouse. Detects a band at 59 kDa. Ideal for mTOR signaling and cell growth studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-S371 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-S371
Species	Human
GeneID (Human)	6198
GeneID	6198
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S371 of human P70S6K1 (NP_001258989.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag