


Anti-Phospho-P90Rsk/Rsk1/Rps6Ka1 (S380) Rabbit Monoclonal Antibody For Wb, Elisa - Q15418

Item Number: CM0001453



Introduction

Phospho-specific rabbit monoclonal antibody targeting RSK1/RPS6KA1 at Ser380. Validated for WB and ELISA in human, mouse, rat. Ideal for studying MAPK signaling, cell proliferation, and survival. Uniprot: Q15418.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk/RSK1/RPS6KA1-S380 Rabbit mAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; Phospho-p90Rsk/RSK1/RPS6KA1-S380
Species	Human
GeneID	6195
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S380 of human p90Rsk/RSK1/RPS6KA1 (Q15418)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag