

Anti-Phospho-Pak2 (S20) Rabbit Monoclonal Antibody For Wb, Elisa - Q13177

Item Number: CM0001110



Introduction

Phospho-specific rabbit monoclonal antibody targeting PAK2 at serine 20. Validated for Western blot and ELISA applications, with cross-reactivity to mouse and rat. Ideal for studying PAK2 signaling in apoptosis and cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PAK2-S20 Rabbit mAb
Remarks/Alias	KNO2; PAK65; PAKgamma; Phospho-PAK2-S20
Species	Human
Gene ID	5062
Immunogen	A synthetic phosphorylated peptide around S20 of human PAK2 (NP_002568.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13177
Protein Weight	58kDa
Shipping	Ice bag