


Anti-Phospho-Paxillin (Y118) Rabbit Monoclonal Antibody For Wb, Elisa - P49023

Item Number: CM0001817



Introduction

Rabbit monoclonal antibody specific for phospho-Paxillin (Tyr118). Validated for WB and ELISA in human, mouse, and rat samples. Ideal for focal adhesion and cell migration studies.

[Learn More](#)

Parameter	Description
Product Name	Phospho-Paxillin-Y118 Rabbit mAb
Remarks/Alias	PXN; paxillin; Phospho-Paxillin-Y118
Species	Human
Gene ID (Human)	5829
Immunogen	A synthetic phosphorylated peptide around Y118 of human Paxillin (P49023)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49023
Protein Weight	65kDa
Shipping	Ice bag