

Anti-Phospho-Pdgfr Alpha-Y754 Rabbit Polyclonal Antibody For Wb, Elisa - P16234

Item Number: CM0001694



Introduction

Phospho-specific rabbit polyclonal antibody targeting PDGFR alpha at Tyr754. Validated for Western blot and ELISA. Detects endogenous levels of PDGFR alpha only when phosphorylated at Tyr754. Human reactive. Useful for studying PDGF signaling and oncology. Predicted molecular weight 123 kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDGFR alpha-Y754 Rabbit pAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2; Phospho-PDGFR alpha-Y754
Species	Human
Gene ID	5156
Immunogen	A synthetic phosphorylated peptide around Y754 of human PDGFRA (NP_006197.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag