


Anti-Phospho-Pdgfrb (Tyr751) Rabbit Polyclonal Antibody For Wb, Elisa - P09619

Item Number: CM0001646



Introduction

Rabbit polyclonal antibody specific to phospho-PDGFRB (Tyr751). Validated for WB and ELISA, reacts with human and mouse. Ideal for studying PDGFRB signaling in cardiovascular development and cell proliferation.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-PDGFRB-Y751 Rabbit pAb
Remarks/Alias	IMF1; KOGS; IBGC4; JTK12; PDGFR; PENTT; CD140B; PDGFR1; PDGFR-1; Phospho-PDGFRB-Y751
Species	Human
GeneID (human)	5159
GeneID	5159
Immunogen	A synthetic phosphorylated peptide around Y751 of human PDGFRB (NP_002600.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09619
Protein Weight	124kDa
Shipping	Ice bag