

Anti-Phospho-Pdha1-S293 Rabbit Polyclonal Antibody For Wb And Elisa - P08559

Item Number: CM0000892



Introduction

High-specificity rabbit polyclonal antibody targeting phosphorylated PDHA1 (serine 293). Validated for Western blot (WB) and ELISA in human, mouse, and rat samples. Ideal for mitochondrial metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDHA1-S293 Rabbit pAb
Remarks/Alias	PDHA; PDHAD; PHE1A; E1alpha; PDHCE1A; Phospho-PDHA1-S293
Species	Human
Gene ID (Human)	5160
Gene ID	5160
Immunogen	Synthetic Peptide / A synthetic phosphorylated peptide around S293 of human PDHA1 (NP_000275.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08559
Protein Weight	43kDa
Shipping	Ice bag