

Anti-Phospho-Perk (T982) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nzj5

Item Number: CM0001768



Introduction

Rabbit monoclonal antibody recognizing phospho-PERK (T982) in human, mouse, rat. Validated for WB, IHC-P, ELISA. Ideal for studying ER stress and unfolded protein response (UPR) signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PERK-T982 Rabbit mAb
Remarks/Alias	PEK; WRS; PERK; Phospho-PERK-T982
Species	Human
GeneID	9451
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T982 of human PERK (NP_004827.4)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ5
Protein Weight	125kDa
Shipping	Ice bag