

Anti-Phospho-Pi3 Kinase P55-Gamma (Pik3R3) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q92569

Item Number: CM0001251

Introduction

Rabbit monoclonal antibody specific for phosphorylated PI3 Kinase p55-gamma (PIK3R3) at Y199. Validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Ideal for insulin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PI3 Kinase p55-Y199 Rabbit mAb
Remarks/Alias	PIK3R3; Phospho-PI3KP55y-Y199
Species	Human
GeneID (Human)	8503
GeneID	8503
Immunogen	A synthetic phosphorylated peptide around Y199 of human PI3 Kinase p55 (NP_001107644.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q92569
Protein Weight	55kDa
Shipping	Ice bag