


Anti-Phospho-Pka C-Alpha (Prkaca)-S339 Polyclonal Antibody For Wb, Ip, Elisa - P17612

Item Number: CM0000437



Introduction

Rabbit polyclonal antibody against phosphorylated PKA C-alpha at Ser339, validated for WB, IP, ELISA in human samples. Ideal for studying cAMP-dependent signaling, cell cycle, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKA C-alpha (PRKACA)-S339 Rabbit pAb
Remarks/Alias	CAFD1; PKACA; PPNAD4; Phospho-PKA C-alpha (PRKACA)-S339
Species	Human
Gene ID (Human)	5566
Gene ID	5566
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S339 of human PKA C-alpha (PRKACA) (NP_002721.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P17612
Protein Weight	41kDa
Shipping	Ice bag