


Anti-Phospho-Pkc Alpha (S319) Rabbit Polyclonal Antibody For Wb And Elisa - P17252

Item Number: CM0000828



Introduction

Phospho-specific rabbit polyclonal antibody targeting PKC alpha phosphorylated at S319, for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for signaling, cancer, and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-S319 Rabbit pAb
Remarks/Alias	AAG6, PKCA, PRKACA, PKCI+/-, PKCalpha, PKC-alpha, Phospho-PKC alpha-S319
Species	Human
Gene ID (Human)	5578
Immunogen	A synthetic phosphorylated peptide around S319 of human PKC alpha (NP_002728.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag