


Anti-Phospho-Pkc Alpha-T638 Rabbit Mab For Wb, Elisa - P17252

Item Number: CM0001410



Introduction

Phospho-specific rabbit monoclonal antibody recognizing human PKC α phosphorylated at Thr638. Ideal for Western blot and ELISA, enabling precise analysis of PKC α activation in cancer, apoptosis, and cardiovascular signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-T638 Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; Phospho-PKC alpha-T638
Species	Human
Gene ID (Human)	5578
Gene ID	5578
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T638 of human PKC alpha.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag