


Anti-Phospho-Pkc Delta-S645 Rabbit Polyclonal Antibody For Wb, Elisa - Q05655

Item Number: CM0001762



Introduction

Phospho-specific rabbit polyclonal antibody targeting PKCδ phosphorylated at S645, validated for WB and ELISA in human, mouse, rat samples. Ideal for apoptosis and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC delta-S645 Rabbit pAb
Remarks/Alias	MAY1; PKCD; ALP53; CVID9; nPKC-delta; Phospho-PKC delta-S645
Species	Human
Gene ID	5580 (Human)
Immunogen	A synthetic phosphorylated peptide around S645 of human PKC delta (NP_006245.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05655
Protein Weight	78kDa
Shipping	Ice bag