


Anti-Phospho-Pkc Delta (T505) Rabbit Polyclonal Antibody For Wb, Elisa - Q05655

Item Number: CM0001442



Introduction

Rabbit polyclonal antibody against human Phospho-PKC delta-T505, validated for WB and ELISA. Targets phosphorylation site T505 of PKC delta (PRKCD, nPKC-delta), a key kinase in apoptosis, cell cycle, and cancer signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC delta-T505 Rabbit pAb
Remarks/Alias	MAY1; PKCD; ALPS3; CVID9; nPKC-delta; Phospho-PKC delta-T505
Species	Human
GeneID (Human)	5580
GeneID	5580
Immunogen	A synthetic phosphorylated peptide around T505 of human PKC delta (NP_006245.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q05655
Protein Weight	78kDa
Shipping	Ice bag