

Anti-Phospho-Pkc Mu (S744/748) Rabbit Polyclonal Antibody For Wb And Elisa - Q15139

Item Number: CM0001245

Introduction

Rabbit polyclonal antibody specific for phosphorylated PKC mu/PRKD1 at Ser744/748, suitable for Western blot and ELISA in human, mouse, and rat samples. Useful for studying PKC mu signaling and phosphorylation events.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC mu-S744/748 Rabbit pAb
Remarks/Alias	PKD; PKCM; CHDED; PRKCM; PKC-MU; Phospho-PKC mu-S744/748
Species	Human
GeneID (Human)	5587
GeneID	5587
Immunogen	Synthetic peptide around S744 & S748 of mouse PKC mu (NP_032884.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15139
Protein Weight	102kDa
Shipping	Ice bag