

Anti-Phospho-Pkc Mu-S910 Rabbit Polyclonal Antibody For Wb, Elisa - Q15139

Item Number: CM0001684

Introduction

Rabbit polyclonal antibody against phosphorylated PKC mu (S910), validated for WB and ELISA in human cells. Suitable for investigating PKD/PRKD1-mediated signaling, Golgi function, and related pathologies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC mu-S910 Rabbit pAb
Remarks/Alias	PKD; PKCM; CHDED; PRKCM; PKC-MU; Phospho-PKC mu-S910
Species	Human
Gene ID	5587
Immunogen	A synthetic phosphorylated peptide around S910 of human PKC mu (NP_002733.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15139
Protein Weight	102kDa
Shipping	Ice bag