


Anti-Phospho-Pkc (Pan) (Gamma)-T514 Rabbit Pab For Wb And Elisa - P17252 / P05771 / Q05655 / Q02156 / P05129 / P24723 / Q04759

Item Number: CM0001019



Introduction

Rabbit polyclonal antibody specific for phosphorylated T514 of PKC pan isoforms (alpha, beta, delta, epsilon, gamma, eta, theta). Validated for WB and ELISA with human, mouse, and rat samples. Ideal for cell signaling, cancer, and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC (pan) (gamma)-T514 Rabbit pAb
Remarks/Alias	PRKCA/PRKCB/PRKCD/PRKCE/PRKCG/PRKCH/PRKCQ; Phospho-PKC (pan) (gamma)-T514
Species	Human
GenelD	5578/5579/5580/5581/5582/5583/5588
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252/P05771/Q05655/Q02156/P05129/P24723/Q04759
Protein Weight	59-67kDa/74-83kDa
Shipping	Ice bag