


Anti-Phospho-Pkc (Pan) (Zeta)-T410 Rabbit Pab - P17252 / P05771 / Q05655 / Q02156 / P05129 / P24723 / P41743 / Q04759 / Q05513

Item Number: CM0000998



Introduction

Phospho-specific rabbit polyclonal antibody targeting T410-phosphorylated PKC zeta and pan-PKC isoforms; validated for WB and ELISA in human, mouse, and rat samples. Suitable for signaling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC (pan) (zeta)-T410 Rabbit pAb
Remarks/Alias	PRKCA/PRKCB/PRKCD/PRKCE/PRKCG/PRKCH/PRKCI/PRKCQ/PRKCZ; Phospho-PKC (pan) (zeta)-T410
Species	Human
GeneID (Human)	5578/5579/5580/5581/5582/5583/5584/5588/5590
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T410 of human PRKCZ (NP_002735.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252/P05771/Q05655/Q02156/P05129/P24723/P41743/Q04759/Q05513
Protein Weight	46kDa/56-68kDa/74-83kDa
Shipping	Ice bag