

Anti-Phospho-Pkc Zeta-T560 Rabbit Monoclonal Antibody For Wb, Elisa - Q05513

Item Number: CM0000699

Introduction

Phospho-specific rabbit monoclonal antibody targeting PKC zeta at Thr560, validated for WB and ELISA. Detects human and cross-reacts with rat. Ideal for studying PKC zeta signaling in cell polarity, NF-kB activation, and insulin response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC zeta-T560 Rabbit mAb
Remarks/Alias	PKC2; PKC-ZETA; Phospho-PKC zeta-T560
Species	Human
GeneID (Human)	5590
GeneID	5590
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T560 of human PKC zeta (Q05513)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q05513
Protein Weight	68 kDa
Shipping	Ice bag