

Anti-Phospho-Pkr/Eif2Ak2-T446 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P19525

Item Number: CM0000797



Introduction

Rabbit polyclonal antibody against phosphorylated PKR/EIF2AK2 (Thr446) for WB, IF/ICC, ELISA. Detects human/mouse target. Ideal for innate immunity, antiviral response, and kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKR/EIF2AK2-T446 Rabbit pAb
Remarks/Alias	PKR; PRKR; DYT33; LEUDEN; EIF2AK1; PPP1R83; Phospho-PKR/EIF2AK2-T446
Species	Human
GeneID (Human)	5610
GeneID	5610
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19525
Protein Weight	62kDa
Shipping	Ice bag