


Anti-Phospho-Pla2G4A-S505 Rabbit Polyclonal Antibody For Wb, Elisa - P47712

Item Number: CM0001127



Introduction

Phospho-specific rabbit polyclonal antibody targeting Ser505 of human cPLA2-alpha (PLA2G4A). Reacts with human, mouse, rat. Validated for WB and ELISA. Ideal for inflammatory pathway research and lipid signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLA2G4A-S505 Rabbit pAb
Remarks/Alias	GURDP; cPLA2; PLA2G4; cPLA2-alpha; Phospho-PLA2G4A-S505
Species	Human
Gene ID (Human)	5321
Gene ID	5321
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47712
Protein Weight	85kDa
Shipping	Ice bag