


Anti-Phospho-Plc Gamma 1 (Plcg1)-Y783 Rabbit Monoclonal Antibody For Wb - P19174

Item Number: CM0001360



Introduction

Rabbit monoclonal antibody recognizing phosphorylated PLC gamma 1 at Tyr783. Validated for WB and ELISA with mouse cross-reactivity. Ideal for RTK signaling and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 1 (PLCG1)-Y783 Rabbit mAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1; Phospho-PLC gamma 1 (PLCG1)-Y783
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	A synthetic phosphorylated peptide around Y783 of human Phospho-PLC gamma 1 (PLCG1)-Y783 (NP_877963.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag