


Anti-Phospho-Plc Gamma 2 (Plcg2)-Y759 Monoclonal Antibody For Wb, If/Icc, Elisa - P16885

Item Number: CM0001409



Introduction

Phospho-specific recombinant rabbit monoclonal antibody against human PLC gamma 2 (PLCG2) phosphorylated at Y759. Validated for WB, IF/ICC, and ELISA; suitable for investigating PLCG2-mediated signal transduction in membrane rafts and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 2 (PLCG2)-Y759 Rabbit mAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2; Phospho-PLC gamma 2 (PLCG2)-Y759
Species	Human
GeneID	5336
Immunogen	Synthetic phosphorylated peptide around Y759 of human PLC gamma 2 (PLCG2) (NP_002652.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag