


Anti-Phospho-Plcg2 (Y1217) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P16885

Item Number: CM0000679



Introduction

Phospho-PLCG2 (Y1217) rabbit monoclonal antibody validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Reacts with human PLC-gamma-2. Ideal for studying tyrosine phosphorylation in B-cell receptor and Fc receptor signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 2 (PLCG2)-Y1217 Rabbit mAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2
Species	Mouse
Gene ID (Human)	5336
Immunogen	Synthetic peptide around Y1217 of human PLC gamma 2 (PLCG2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag