


Anti-Phospho-Plcg2 (Y759) Rabbit Polyclonal Antibody For Wb, Elisa - P16885

Item Number: CM0001459



Introduction

Phospho-specific rabbit polyclonal antibody against human PLCG2 (Y759) for WB and ELISA applications. This antibody targets a crucial enzyme in transmembrane signaling that produces second messengers DAG and IP3. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 2 (PLCG2)-Y759 Rabbit pAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2; Phospho-PLC gamma 2 (PLCG2)-Y759
Species	Human
Gene ID	5336
Immunogen	A synthetic phosphorylated peptide around Y759 of human PLC gamma 2 (PLCG2) (NP_002652.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag