


Anti-Phospho-Plk1-T210 Rabbit Monoclonal Antibody For Wb And Elisa - P53350

Item Number: CM0000858



Introduction

Rabbit monoclonal antibody against human PLK1 phosphorylated at Thr210. Generated using a synthetic phospho-peptide immunogen. Ideal for cell cycle and mitosis studies using Western blot and ELISA.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLK1-T210 Rabbit mAb
Remarks/Alias	PLK; STPK13; Phospho-PLK1-T210
Species	Human
Gene ID	5347 (Human)
Immunogen	A synthetic phosphorylated peptide around T210 of human PLK1 (P53350)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P53350
Protein Weight	68 kDa
Shipping	Ice bag