


Anti-Phospho-Polr2A Ctd-S2 Polyclonal Antibody For Wb, Ip, Elisa - P24928

Item Number: CM0000948



Introduction

Rabbit polyclonal antibody recognizing human POLR2A phosphorylated at Ser2 of the CTD repeat (YSPTSPS). Validated for Western blot, immunoprecipitation, and ELISA. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-S2 Rabbit pAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplILS; hsRPB1; hRPB220; Phospho-POLR2A CTD-S2
Species	Human
Gene ID (Human)	5430
Gene ID	5430
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S2 of human POLR2A CTD repeat YSPTSPS.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag