


Anti-Phospho-Polr2A Ctd-S5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24928

Item Number: CM0001299



Introduction

High-quality anti-phospho-POLR2A CTD-S5 rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets RNA polymerase II CTD Ser5 phosphorylation, essential for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-S5 Rabbit pAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplLS; hsRPB1; hRPB220; Phospho-POLR2A CTD-S5
Species	Human
Gene ID (Human)	5430
Gene ID	5430
Immunogen	Synthetic Peptide: A phospho specific peptide corresponding to residues surrounding S5 of human POLR2A CTD repeat YSPTSPS.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag