

Anti-Phospho-Ppp1Ca (T320) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P62136

Item Number: CM0000476



Introduction

Rabbit polyclonal antibody against phosphorylated PPP1CA (Thr320) validated in WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying PP1 α catalytic subunit phosphorylation and phosphatase signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PPP1CA-T320 Rabbit pAb
Remarks/Alias	PP1A; PP-1A; PPP1A; PP1alpha; Phospho-PPP1CA-T320
Species	Human
GeneID (Human)	5499
GeneID	5499
Immunogen	A synthetic phosphorylated peptide around T320 of human PPP1CA (NP_002699.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62136
Protein Weight	38kDa
Shipping	Ice bag