


Anti-Phospho-Ppp1R12A (S507) Rabbit Polyclonal Antibody For Wb, Elisa - O14974

Item Number: CM0001861



Introduction

Phospho-specific rabbit polyclonal antibody targeting Ser507 of human PPP1R12A (MYPT1), validated for Western blot and ELISA. Ideal for studying myosin phosphatase regulation and smooth muscle signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PPP1R12A-S507 Rabbit pAb
Remarks/Alias	MBS; GUBS; M130; MYPT1; Phospho-PPP1R12A-S507
Species	Human
GeneID (human)	4659
GeneID	4659
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14974
Protein Weight	115kDa
Shipping	Ice bag