

Anti-Phospho-Pten-S380/T382/383 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P60484

Item Number: CM0001623



Introduction

Phospho-PTEN-S380/T382/383 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat PTEN; targets phosphorylated tumor suppressor PTEN in PI3K-AKT signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PTEN-S380/T382/383 Rabbit pAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; Phospho-PTEN-S380/T382/383
Species	Human
Gene ID (Human)	5728
Gene ID	5728
Immunogen	Synthetic phosphorylated peptide around S380, T382, T383 of human PTEN (NP_000305.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag