

Anti-Phospho-Pten-S380/T382/T383 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P60484

Item Number: CM0001402



Introduction

Phospho-PTEN-S380/T382/T383 rabbit monoclonal antibody detects endogenous PTEN phosphorylated at S380/T382/T383. Validated for WB, IF-P, IHC-P, ELISA with human, mouse, rat reactivity. Ideal for PI3K-AKT pathway and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PTEN-S380/T382/T383 Rabbit mAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; Phospho-PTEN-S380/T382/T383
Species	Human
GeneID (human)	5728
GeneID	5728
Immunogen	A synthetic phosphorylated peptide around S380 & T383 of human Phospho-PTEN-S380/T382/T383 (NP_000305.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag