

Anti-Phospho-Pten (Ser380) Rabbit Polyclonal Antibody For Wb And Elisa - P60484

Item Number: CM0000651



Introduction

Rabbit polyclonal antibody raised against a synthetic phosphopeptide around S380 of human PTEN. Validated for Western blot (WB) and ELISA. Detects human, mouse, and rat PTEN phosphorylated at S380. Ideal for studying PTEN phosphorylation in PI3K/AKT signaling and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PTEN-S380 Rabbit pAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; Phospho-PTEN-S380
Species	Human
Gene ID (Human)	5728
Gene ID	5728
Immunogen	A synthetic phosphorylated peptide around S380 of human PTEN (NP_000305.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag