


Anti-Phospho-Ptk2B (Y402) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14289

Item Number: CM0001514



Introduction

Rabbit polyclonal antibody specific for PTK2B phosphorylated at tyrosine 402. Validated for Western blot, IF/ICC, and ELISA. Detects PTK2B in human, mouse, and rat samples. Suitable for studying integrin signaling, cell migration, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PTK2B-Y402 Rabbit pAb
Remarks/Alias	PKB; PTK; CAKB; FAK2; PYK2; CADTK; FADK2; RAFTK; Phospho-PTK2B-Y402
Species	Human
GeneID (Human)	2185
GeneID	2185
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y402 of human PTK2B (NP_004094.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14289
Protein Weight	116kDa
Shipping	Ice bag