

Anti-Phospho-Ptptra-Y789 Polyclonal Antibody For Wb, Elisa - P18433

Item Number: CM0000507



Introduction

High-quality anti-Phospho-PTPRA-Y789 rabbit polyclonal antibody for WB and ELISA. Cross-reacts with human, mouse, rat PTPRA. Useful in integrin-mediated focal adhesion and cell junction research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-PTPRA-Y789 Rabbit pAb
Remarks/Alias	LRP; HLPRA; PTPA; HEPTP; HPTPA; RPTPA; PTPRL2; HPTPalpha; R-PTP-alpha; Phospho-PTPRA-Y789
Species	Human
GeneID (Human)	5786
GeneID	5786
Immunogen	Synthetic phosphorylated peptide around Y789 of human PTPRA (NP_543030.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18433
Protein Weight	91kDa
Shipping	Ice bag