

Anti-Phospho-Raf1-S259 Rabbit Monoclonal Antibody For Wb And Elisa - P04049

Item Number: CM0001530



Introduction

Phospho-Raf1-S259 rabbit monoclonal antibody for Western blot (WB) and ELISA; detects human and rat RAF1 phosphorylated at serine 259. Suitable for studying MAPK/ERK signaling, cell proliferation, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raf1-S259 Rabbit mAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Phospho-Raf1-S259
Species	Human
GeneID (Human)	5894
GeneID	5894
Immunogen	A synthetic phosphorylated peptide around S259 of human RAF1 (P04049)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag