

Anti-Phospho-Raf1-S621 Rabbit Monoclonal Antibody For Wb, Elisa - P04049

Item Number: CM0000569



Introduction

Phospho-specific rabbit monoclonal antibody targeting human Raf1 (c-Raf) phosphorylated at serine 621. Validated for WB and ELISA; cross-reacts with rat. Ideal for MAPK/ERK signaling studies and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raf1-S621 Rabbit mAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Phospho-Raf1-S621
Species	Human
GeneID (Human)	5894
GeneID	5894
Immunogen	A synthetic phosphorylated peptide around S621 of human Phospho-Raf1-S621 (P04049)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag