

Anti-Phospho-Raf1 (Ser621) Rabbit Polyclonal Antibody - P04049

Item Number: CM0001876



Introduction

Rabbit polyclonal antibody specific for phosphorylated Raf1 at Ser621. Validated for Western blot, immunoprecipitation, and ELISA. Recognizes human Raf1, a key kinase in MAPK signaling and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raf1-S621 Rabbit pAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Phospho-Raf1-S621
Species	Human
GeneID (human)	5894
GeneID	5894
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S621 of human Phospho-Raf1-S621 (NP_002871.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag