

Anti-Phospho-Raptor-Ser792 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N122

Item Number: CM0000850



Introduction

Rabbit polyclonal antibody targeting phosphorylated Raptor (Ser792) for Western blot and ELISA. Validated in human and mouse, this antibody monitors mTORC1 pathway activation; suitable for research in cancer metabolism and cell growth.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raptor-Ser792 Rabbit pAb
Remarks/Alias	KOG1; Mip1; Phospho-Raptor-Ser792
Species	Human
GeneID (Human)	57521
GeneID	57521
Immunogen	A synthetic phosphorylated peptide around S792 of human RPTOR (NP_065812.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N122
Protein Weight	149kDa
Shipping	Ice bag