


Anti-Phospho-Rb (S795) Rabbit Polyclonal Antibody For Wb - P06400

Item Number: CM0001548



Introduction

Phospho-RB-S795 rabbit polyclonal antibody specific for retinoblastoma protein phosphorylated at Ser795, validated for WB, IHC-P, ELISA, reacts with human, mouse, rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S795 Rabbit pAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S795
Species	Human
GeneID (Human)	5925
GeneID	5925
Immunogen	A synthetic phosphorylated peptide around S795 of human RB (NP_000312.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag