


Anti-Phospho-Rb-S807/811 Rabbit Monoclonal Antibody For Wb - P06400

Item Number: CM0001031



Introduction

Phospho-RB-S807/811 Rabbit monoclonal antibody detects phosphorylated RB at S807/811 in human and mouse, suitable for WB, IP, IF/ICC, IF-P, and ELISA, supporting cell cycle and tumor suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S807/811 Rabbit mAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S807/811
Species	Human
Gene ID (Human)	5925
Gene ID	5925
Immunogen	A synthetic phosphorylated peptide around S807 & S811 of human RB (NP_000312.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag