


Anti-Phospho-Rb1 (S780) Rabbit Monoclonal Antibody For Wb, Elisa - P06400

Item Number: CM0001744



Introduction

Phospho-specific rabbit monoclonal antibody against RB1 (S780). Validated for Western blot and ELISA with human and rat samples. Supports investigation of cell cycle control and tumor suppressor pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S780 Rabbit mAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S780
Species	Human
Gene ID (Human)	5925
Gene ID	5925
Immunogen	Synthetic peptide A phospho synthetic peptide corresponding to residues surrounding S780 of Human RB.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag