


Anti-Phospho-Rb1 (Ser807/Ser811) Rabbit Polyclonal Antibody For Wb, Elisa - P06400

Item Number: CM0001225



Introduction

Rabbit polyclonal antibody targeting phosphorylated RB1 at Ser807/Ser811. Suitable for Western blot and ELISA with cross-reactivity to human, mouse, and rat. Key regulator of G1/S cell cycle transition.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S807/811 Rabbit pAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S807/811
Species	Human
GeneID (human)	5925
GeneID	5925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S807 & S811 of human RB (NP_000312.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag