


Anti-Phospho-Ripk1 (S166) Rabbit Monoclonal Antibody For Wb, Elisa - Q13546

Item Number: CM0001043



Introduction

Anti-Phospho-RIPK1 (S166) rabbit monoclonal antibody validated for WB and ELISA. Detects human RIPK1 phosphorylated at serine 166. Ideal for cell death and inflammatory signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RIPK1-S166 Rabbit mAb
Remarks/Alias	RIP; RIP1; AIEFL; IMD57; RIP-1; Phospho-RIPK1-S166
Species	Human
Gene ID (Human)	8737
Gene ID	8737
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S166 of human RIPK1 (NP_003795.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13546
Protein Weight	75kDa
Shipping	Ice bag