


## Anti-Phospho-Ripk1 (S166) Rabbit Polyclonal Antibody For Wb And Elisa - Q13546

Item Number: CM0001155



### Introduction

Rabbit polyclonal antibody recognizing phosphorylated RIPK1 at Ser166 (human, mouse). Validated for Western blot and ELISA. Ideal for studying TNF-mediated apoptosis, necroptosis, and inflammatory signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                           |
|------------------|---------------------------------------------------------------------------------|
| Product Name     | Phospho-RIPK1/RIP-S166 Rabbit pAb                                               |
| Remarks/Alias    | RIP; RIP1; AIEFL; IMD57; RIP-1; Phospho-RIPK1/RIP-S166                          |
| Species          | Human                                                                           |
| Gene ID          | 8737                                                                            |
| Immunogen        | A synthetic phosphorylated peptide around S166 of human RIPK1/RIP (NP_003795.2) |
| Source           | Rabbit                                                                          |
| Category         | Phosphorylated Antibodies                                                       |
| Application      | WB, ELISA                                                                       |
| Cross Reactivity | Human, Mouse                                                                    |
| SWISS            | Q13546                                                                          |
| Protein Weight   | 76 kDa                                                                          |
| Shipping         | Ice bag                                                                         |