


Anti-Phospho-Rpa2-Ser33 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15927

Item Number: CM0000857



Introduction

Phospho-specific rabbit monoclonal antibody against human RPA2 phosphorylated at Ser33, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat RPA2. Ideal for studying DNA replication, repair, and damage signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RPA2-Ser33 Rabbit mAb
Remarks/Alias	REPA2; RPA32; RP-A p32; RP-A p34
Species	Human
GeneID	6118
Immunogen	Synthetic phosphorylated peptide around Ser33 of human RPA2 (NP_002937.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P15927
Protein Weight	29 kDa
Shipping	Ice bag